


Product Datasheet


C18orf25 antibody (orb325690)

Description	Rabbit polyclonal antibody to C18orf25
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human C18orf25
Target	C18orf25
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: STSSDDDEEVSGSKTITAEIPGHLDPGFLASDKTSAGNAPLNNEINIA
Clonality	Polyclonal
MW	37kDa
Uniprot ID	Q96B23
NCBI	NP_001008240
Expiration Date	12 months from date of receipt.

C18orf25



Western blot analysis of human Fetal Liv...



Western blot analysis of 721_B cell lysa...

C18orf25



Western blot analysis of human Fetal Hea...