

Product Datasheet

DnaI4 antibody (orb325718)

Description	Rabbit polyclonal antibody to DnaI4
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	DnaI4
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ESAAKMIKETMDKKFGSSWHVVIGEGFGFEITHEVKNLLYLYFGGTLAVC
Clonality	Polyclonal
MW	12kDa
Uniprot ID	Q9DCM4
NCBI	NP_059498
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Small Intestine (Small Intestine) showing a single band at approximately 12 kDa, corresponding to the DnaI4 protein.