


Product Datasheet

CCDC144NL antibody (orb325908)

Description	Rabbit polyclonal antibody to CCDC144NL
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human MGC87631
Target	CCDC144NL
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VDKRD RKKSIIQLVPEYKEKQTPE SL PQNNNPAAPSQAEGGEGGVACGTV
Clonality	Polyclonal
MW	24
Uniprot ID	Q6NUI1
NCBI	NP_001004306
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...