

Product Datasheet

GRHL2 antibody (FITC) (orb32606)

Description	GRHL2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Grainyhead-like protein 2 homolog protein (1-265AA)
Target	GRHL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6ISB3
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.