

Product Datasheet

NOTCH2NL antibody (orb326061)

Description	Rabbit polyclonal antibody to NOTCH2NL
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Target	NOTCH2NLA
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LYVPCAPSPCVNGGTCRQTGDFTFECNCLPETVRRGTELWERDREVVWNGK
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q7Z3S9
NCBI	NP_982283
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat Whole Cell lysate.