

Product Datasheet

ERVFRD-1 antibody (orb326066)

Description	Rabbit polyclonal antibody to ERVFRD-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human ERVFRD-1
Target	ERVFRD-1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PLVSLLLLLLFGPCLLNLTQFVSSRLQAIKLQTNLSAGRHPNRIQESPF
Clonality	Polyclonal
MW	59kDa
Uniprot ID	P60508
NCBI	NP_997465
Expiration Date	12 months from date of receipt.



Western blot analysis of human Thymus Tu...