

Product Datasheet

KIAA1841 antibody (orb32609)

Description	KIAA1841 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Uncharacterized protein KIAA1841 protein (401-707AA)
Target	KIAA1841
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6NSI8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.