


Product Datasheet

MUC3B antibody (orb326253)

Description	Rabbit polyclonal antibody to MUC3B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human MUC3B
Target	MUC3B
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KTTLKEGLQNASQDANSCQDSQTLCKPDSIKVNNNSKTELTPEAICRRA
Clonality	Polyclonal
MW	107kDa
Uniprot ID	Q9H195
NCBI	XP_001125753
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7 cell lysat...