



Product Datasheet

ERLIN2 antibody (orb326257)


Description	Rabbit polyclonal antibody to ERLIN2
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ERLIN2
Target	ERLIN2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: ASNSKIYFGKDIPNMFMDMSAGSVSKQFEGGLADKLSFGLEDEPLETATKEN
Clonality	Polyclonal
MW	38kDa
Uniprot ID	O94905
NCBI	NP_009106
Expiration Date	12 months from date of receipt.



Western blot analysis of human Fetal Lung...



Western blot analysis of human Fetal Heart...



Western blot analysis of human Fetal Brain...