


Product Datasheet

UBQLN2 antibody (orb326289)

Description	Rabbit polyclonal antibody to UBQLN2
Species/Host	Rabbit
Reactivity	Canine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human UBQLN2
Target	UBQLN2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LKDQDTLIQHGIHDGLTVHLVIKSQNRPPQGQSTQPSNAAGTNTTSASTPR
Clonality	Polyclonal
MW	66kDa
Uniprot ID	Q9UHD9
NCBI	NP_038472
Expiration Date	12 months from date of receipt.

UBQLN2



Western blot analysis of human HeLa tiss...



Western blot analysis of HepG2 cell lysa...