


Product Datasheet

NAA38 antibody (orb326340)

Description	Rabbit polyclonal antibody to NAA38
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen for Anti-NAA38 antibody is: synthetic peptide directed towards the C-terminal region of Human LSMD1
Target	NAA38
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9BRA0
Expiration Date	12 months from date of receipt.



Western blot analysis of Uterus Tumor ti...