

Product Datasheet

ARHGAP11A antibody (Biotin) (orb32642)

Description	ARHGAP11A antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Rho GTPase-activating protein 11A protein (1-300AA)
Target	ARHGAP11A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6P4F7
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.