

Product Datasheet

PPP6R3 antibody (orb326422)

Description	Rabbit polyclonal antibody to PPP6R3
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	PPP6R3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PSQEEDRHSNASQSLCEIVRLSRDQMLQIQNSTEPDPLLATLEKQEIEQ
Clonality	Polyclonal
MW	92kDa
Uniprot ID	Q5H9R7
NCBI	NP_060782
Expiration Date	12 months from date of receipt.



Western blot analysis of 721_B Whole Cel...