

Product Datasheet

CRISPLD1 antibody (orb326444)

Description	Rabbit polyclonal antibody to CRISPLD1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	WB
Target	CRISPLD1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9H336
Expiration Date	12 months from date of receipt.



Western blot analysis of Fetal Liver tis...