

Product Datasheet

SAA2 antibody (orb326476)

Description	Rabbit polyclonal antibody to SAA2
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Porcine, Rat, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Target	SAA2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GAWAAEVISNARENIQRLTGRGAEDSLADQAANKWGRSGRDPNHFRPAGL
Clonality	Polyclonal
MW	14 kDa
Uniprot ID	P0DJ19
NCBI	NP_110381
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T Whole Cell...