

Product Datasheet

SYNGAP1 antibody (orb326496)

Description	Rabbit polyclonal antibody to SYNGAP1
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Target	SYNGAP1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: EYKLKEYSKSMDESRLDRVKEYEEIHSLKERLHMSNRKLEEYERRLLSQ
Clonality	Polyclonal
MW	148kDa
Uniprot ID	Q96PV0
NCBI	NP_006763
Expiration Date	12 months from date of receipt.



Western blot analysis of Fetal kidney tissue...