

Product Datasheet

CC2D2A antibody (orb326553)

Description	Rabbit polyclonal antibody to CC2D2A
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human CC2D2A
Target	CC2D2A
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	E7EP21
Expiration Date	12 months from date of receipt.



Western blot analysis of MDA-MB-435S Who...