


Product Datasheet

PKHB1 antibody (orb326560)

Description	Rabbit polyclonal antibody to PKHB1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human PKHB1
Target	PLEKHB1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	Batch dependent within range: 100 ul at 0.5 - 1 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LWLDGTLGYYHDETAQDEEDRVLIHFNVVDIKIGPECHDVQPPEGRSRDG
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q9UF11
NCBI	NP_067023
Expiration Date	12 months from date of receipt.



Western blot analysis of human 786-0 cells treated with PKHB1 antibody.