

Product Datasheet

C16orf93 antibody (orb326846)

Description	Rabbit polyclonal antibody to C16orf93
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human C16orf93
Target	CCDC189
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: IKTQVNKELEQLQGLVEERLKASEERLSSKLTALERPFQLPPGKGKSKTK
Clonality	Polyclonal
MW	32kDa
Uniprot ID	A1A4V9
NCBI	NP_001182549
Expiration Date	12 months from date of receipt.

