


Product Datasheet

CCDC42B antibody (orb326921)

Description	Rabbit polyclonal antibody to CCDC42B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human CCDC42B
Target	CFAP73
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: QKERELKGSFIRFDKFLQDSEARRNRALRRAAEERHQAGRREVEALRLWT
Clonality	Polyclonal
MW	33kDa
Uniprot ID	A6NFT4
NCBI	XP_005253936
Expiration Date	12 months from date of receipt.



Western blot analysis of human HeLa cells transfected with CCDC42B. The blot shows a single band at approximately 33 kDa, indicating successful transfection and expression of the target protein.