

Product Datasheet

C6orf99 antibody (orb327152)

Description	Rabbit polyclonal antibody to C6orf99
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen for Anti-C6orf99 antibody is: synthetic peptide directed towards the middle region of Human CF099
Target	C6orf99
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: PGSFFLCKIRECVLNRYRFQLQHPGFQHYLQSSGRRDRGRSEDKKPLEAGV
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	Q4VX62
Expiration Date	12 months from date of receipt.



Western blot analysis of human HCT116 cells treated with...
Wh...