


Product Datasheet

NSMF antibody (orb327366)

Description	Rabbit polyclonal antibody to NSMF
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human NSMF
Target	NSMF
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: SQSHPENRNGADHLLADAYSGHDGSPQMPPAPQNKRRSLVSNCGYEGSL
Clonality	Polyclonal
MW	58kDa
Uniprot ID	Q6X4W1
Expiration Date	12 months from date of receipt.



Western blot analysis of human 293T cells. Whole cell lysates were probed with NSMF antibody (orb327366).