

Product Datasheet

EAF6 antibody (orb327393)

Description	Rabbit polyclonal antibody to EAF6
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N-terminal region of Human EAF6
Target	MEAF6
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: GWDRLTNQKNSNSKNDRRNRKFKEAERLFSKSSVTSAAAVSALAGVQDQ
Clonality	Polyclonal
MW	21kDa
Uniprot ID	Q9HAF1
Expiration Date	12 months from date of receipt.



Western blot analysis of human U937 cells. Whole cell lysates were probed with EAF6 antibody.