



Product Datasheet

PDL1 antibody (orb32884)

Description	Rabbit polyclonal antibody to PDL1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	CD274
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	33275
Uniprot ID	Q9NZQ7
NCBI	NP_054862.1 , NP_001254635.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of A549 whole cell...



Western blot analysis of A549 whole cell...