


Product Datasheet

beta 1 Adrenergic Receptor antibody [Out of stock] (orb329571)

Description	Rabbit polyclonal antibody to beta 1 Adrenergic Receptor
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C-terminal region of Human ADRB1
Target	ADRB1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08588
Expiration Date	12 months from date of receipt.



Western blot analysis of human THP-1 tis...