



Product Datasheet

POU2F3 antibody (orb329669)


Description	Rabbit polyclonal antibody to POU2F3
Species/Host	Rabbit
Reactivity	Canine, Guinea pig, Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human POU2F3
Target	POU2F3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: KMSGDVADSTDARSTLSQVEPGNDRKGLDFNRQIKTEDLSDSLQQTLSHR
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	Q9UKI9
NCBI	NP_055167
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse to...



Western blot analysis of HepG2 cell lysa...



Immunohistochemical staining of human Mu...