


Product Datasheet

POU2F3 antibody (orb329670)

Description	Rabbit polyclonal antibody to POU2F3
Species/Host	Rabbit
Reactivity	Canine, Guinea pig, Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human POU2F3
Target	POU2F3
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UKI9
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 Whole Cel...



Immunohistochemical staining of mouse to...