

Product Datasheet

AKAP8L antibody (orb329673)

Description	Rabbit polyclonal antibody to AKAP8L
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human AKAP8L
Target	AKAP8L
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9ULX6
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...