


Product Datasheet

PPP1R13L antibody (orb329694)

Description	Rabbit polyclonal antibody to PPP1R13L
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PPP1R13L
Target	PPP1R13L
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: AQSVPLEEVARVLAEIPRPLKRRGSMEQAPAVALPPTHKKQYQQIISRL
Clonality	Polyclonal
MW	89kDa
Uniprot ID	Q8WUF5
NCBI	NP_006654
Expiration Date	12 months from date of receipt.


PPP1R13L




Green: PPP1R13L
Blue: DAPI

Immunohistochemical staining of human A549 cells

PPP1R13L



Western blot analysis of PPP1R13L antibody



Western blot analysis of Transfected 293T cells