



Product Datasheet

GTF2IRD1 antibody (orb329740)

Description	Rabbit polyclonal antibody to GTF2IRD1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human GTF2IRD1
Target	GTF2IRD1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TFGSQNLERILAVADKIKFTVTRPFQGLIPKPEDDDANRLGEKVILREQV
Clonality	Polyclonal
MW	106kDa
Uniprot ID	Q9UHL9
NCBI	NP_057412
Expiration Date	12 months from date of receipt.



Western blot analysis of Transfected 293...



Immunohistochemical staining of human ki...