

Product Datasheet


PARP2 antibody (orb329744)

Description	Rabbit polyclonal antibody to PARP2
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human PARP2
Target	PARP2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LLDLFEVEKDGEKEAFREDLHNRMLLWHGSRMSNWVGILSHGLRIAHPEA
Clonality	Polyclonal
MW	66 kDa
Uniprot ID	Q9UGN5
NCBI	NP_005475
Expiration Date	12 months from date of receipt.



Western blot analysis of Raji cell lysat...

PARP2



Immunofluorescence analysis of rat thyro...