



Product Datasheet

C4BPB antibody (orb329748)


Description	Rabbit polyclonal antibody to C4BPB
Species/Host	Rabbit
Reactivity	Canine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human C4BPB
Target	C4BPB
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: CCLMVAWRVSASDAEHCPVPPVDNSIFVAKEVEGQILGTYVCIKGYHLV
Clonality	Polyclonal
MW	28kDa
Uniprot ID	P20851
NCBI	NP_000707
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human kidney tissue.



Western blot analysis of human 721_B tissue lysate.



Western blot analysis of human Liver tissue lysate.