



Product Datasheet

RGS13 antibody (orb329758)

Description	Rabbit polyclonal antibody to RGS13
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human RGS13
Target	RGS13
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: WSRISRAKKLYKIYIQQSPREINIDSSTRETIIRNIQEPTETCFEEAQK
Clonality	Polyclonal
MW	19kDa
Uniprot ID	O14921
NCBI	NP_002918
Expiration Date	12 months from date of receipt.



Western blot analysis of Daudi cell lysa...



Immunohistochemical staining of human ki...