


Product Datasheet

EZH1 antibody [Out of stock] (orb329768)

Description	Rabbit polyclonal antibody to EZH1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human EZH1
Target	EZH1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q92800
Expiration Date	12 months from date of receipt.



Western blot analysis of HCT116 tissue u...



Western blot analysis of human brain tis...