



Product Datasheet

KCNMA1 antibody (orb329808)


Description	Rabbit polyclonal antibody to KCNMA1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Sheep, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the c terminal region of human KCNMA1
Target	KCNMA1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: CFGYRLRDAHLSTPSQCTKRYVITNPPYEFELVPTDLIFCLMQFDHNAG
Clonality	Polyclonal
MW	131kDa
Uniprot ID	Q12791
NCBI	NP_001014797
Expiration Date	12 months from date of receipt.



Western blot analysis of human Liver tis...



Immunohistochemical staining of human Pr...



Immunohistochemical staining of human Pi...