


Product Datasheet

Kcnq3 antibody (orb329812)

Description	Rabbit polyclonal antibody to Kcnq3
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide corresponding to a region of Mouse
Target	Kcnq3
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: TPKHKKSQKGSFTYPSQQSPRNEPYVARAATSETEDQSMMGKFVKVERQ
Clonality	Polyclonal
MW	97kDa
Uniprot ID	Q8K3F6
NCBI	NP_690887
Expiration Date	12 months from date of receipt.



Western blot analysis of Kcnq3 antibody



Western blot analysis of mouse Heart tis...