



Product Datasheet

CACNA1G antibody (orb329837)

Description	Rabbit polyclonal antibody to CACNA1G
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human CACNA1G
Target	CACNA1G
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VSRTHSLPNDSYMCRHGSTAEGPLGHRGWGLPKAQSGSVLSVHSQPADTS
Clonality	Polyclonal
MW	241 kDa
Uniprot ID	O43497
NCBI	NP_938202
Expiration Date	12 months from date of receipt.



Western blot analysis of human Spleen ti...



Immunohistochemical staining of human br...