


Product Datasheet

KCTD11 antibody (orb329839)

Description	Rabbit polyclonal antibody to KCTD11
Species/Host	Rabbit
Reactivity	Equine, Goat, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human KCTD11
Target	KCTD11
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: RLGRDLPRGYGETALLRAEADFYQIRPLLDALRELEASESQGTPAPTAALL
Clonality	Polyclonal
MW	26kDa
Uniprot ID	Q693B1
NCBI	NP_001002914
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat cell lys...