


Product Datasheet

PIAS2 antibody (orb329847)


Description	Rabbit polyclonal antibody to PIAS2
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of human PIAS2
Target	PIAS2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	1.0 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FLDSLTSPLTASSTSVTTTSSHESSTHVSSSSSRSETGVITSSGSNIPDI
Clonality	Polyclonal
MW	68kDa
Uniprot ID	O75928
NCBI	NP_004662
Expiration Date	12 months from date of receipt.

PIAS2



Western blot analysis of human HepG2 tis...

PIAS2



Western blot analysis of human 721_B tis...



Western blot analysis of HepG2 cell lysis...