



Product Datasheet

beta 1 Adrenergic Receptor antibody (orb329873)


Description	Rabbit polyclonal antibody to beta 1 Adrenergic Receptor
Species/Host	Rabbit
Reactivity	Canine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ADRB1
Target	ADRB1
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: CTVWAISALVSFLPILMHWWRAESDEARRCYNDPKCCDFVTNRAYAIASS
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P08588
NCBI	NP_000675
Expiration Date	12 months from date of receipt.



Western blot analysis of human brain tis...



Western blot analysis of Wild type mouse...



Western blot analysis of human Placenta ...