


Product Datasheet

Pias2 antibody (orb329894)

Description	Rabbit polyclonal antibody to Pias2
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the C terminal region of mouse Pias2
Target	Pias2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LSLIPVDPQYCPMFLDSLTSPLTASSTSVTTTSPHESSTHVSSSSSRSE
Clonality	Polyclonal
MW	54kDa
Uniprot ID	Q8C5D8
NCBI	AAB96678
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse Pa...



Western blot analysis of human heart tis...



Western blot analysis of NIH/3T3 cell ly...