


Product Datasheet

E2f7 antibody (orb329925)


Description	Rabbit polyclonal antibody to E2f7
Species/Host	Rabbit
Reactivity	Canine, Equine, Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of mouse E2f7
Target	E2f7
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MEVNCLTLKDLISPRQTRLDFAIEDAENAQKENIFVDRSRMTPKTPMKNE
Clonality	Polyclonal
MW	99kDa
Uniprot ID	Q8BSQ3
NCBI	NP_848724
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse st...



Western blot analysis of human SP2/0 Who...



Western blot analysis of SP2/0 cell lysa...