


Product Datasheet

MEFV antibody (orb329947)

Description	Rabbit polyclonal antibody to MEFV
Species/Host	Rabbit
Reactivity	Equine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human MEFV
Target	MEFV
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MAKTPSDHLLSTLEELVPYDFEKFKFLQNTSVQKEHSRIPRSQIQRARP
Clonality	Polyclonal
MW	86kDa
Uniprot ID	O15553
NCBI	NP_000234
Expiration Date	12 months from date of receipt.



Western blot analysis of human Fetal Brain lysate.



Western blot analysis of HepG2 cell lysate.