

Product Datasheet

NR4A2 antibody (orb330014)

Description	Rabbit polyclonal antibody to NR4A2
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human NR4A2
Target	NR4A2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	Batch dependent within range: 100 ul at 0.5 - 1 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: FSCIAALAMVTERHGLKEPKRVEELQNKIVNCLKDHVTFNNGGLNRPNYL
Clonality	Polyclonal
MW	66kDa
Uniprot ID	P43354
NCBI	NP_006177
Expiration Date	12 months from date of receipt.



Western blot analysis of human Liver tis...