



Product Datasheet

PPARGC1A antibody (orb330034)


Description	Rabbit polyclonal antibody to PPARGC1A
Species/Host	Rabbit
Reactivity	Canine, Equine, Goat, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human PPARGC1A
Target	PPARGC1A
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: MAWDMCNQDSESVWSDIECAALVGEDQPLCPDLPELDLSELDVNDLDTDS
Clonality	Polyclonal
MW	91 kDa
Uniprot ID	Q9UBK2
NCBI	NP_037393
Expiration Date	12 months from date of receipt.



Western blot analysis of PPARGC1A antibody...



Western blot analysis of ACHN cell lysate...



Western blot analysis of transfected SH- ...