



Product Datasheet

ZEB2 antibody (orb330041)


Description	Rabbit polyclonal antibody to ZEB2
Species/Host	Rabbit
Reactivity	Equine, Guinea pig, Human, Mouse, Porcine, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	The immunogen is a synthetic peptide directed towards the middle region of human ZEB2
Target	ZEB2
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: LGRQDGDEEFEEEESEENKSMDDTDPETIRDEEETGDHSMDDSSSEDGKME
Clonality	Polyclonal
MW	25 kDa
Uniprot ID	O60315
NCBI	NP_055610
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human Sm...



Western blot analysis of human H9 tissue...



Western blot analysis of HeLa cell lysat...