


Product Datasheet

ZMIZ1 antibody (orb330053)

Description	Rabbit polyclonal antibody to ZMIZ1
Species/Host	Rabbit
Reactivity	Canine, Equine, Guinea pig, Human, Mouse, Rat, Zebrafish
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human ZMIZ1
Target	ZMIZ1
Preservatives	Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8C up to 1 week. For long term storage, store at -20C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9ULJ6
Expiration Date	12 months from date of receipt.



Western blot analysis of ZMIZ1 antibody



Western blot analysis of 293T cell lysat...