



Product Datasheet

VCAM1 antibody (orb33009)

Description	Rabbit polyclonal antibody to VCAM1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	VCAM1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	81276
Uniprot ID	P19320
NCBI	NP_001069.1 , NP_542413.1 , NP_001186763.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell line...