

Product Datasheet

alpha 2 Macroglobulin antibody (orb33013)

Description	Rabbit polyclonal antibody to alpha 2 Macroglobulin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	A2M
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	163291
Uniprot ID	P01023
NCBI	NP_000005.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

