


Product Datasheet

SIGLEC6 antibody (orb330211)

Description	Rabbit polyclonal antibody to SIGLEC6
Species/Host	Rabbit
Reactivity	Equine, Human, Porcine
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	The immunogen is a synthetic peptide directed towards the N terminal region of human SIGLEC6
Target	SIGLEC6
Preservatives	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Form/Appearance	Liquid. Purified antibody supplied in 1x PBS buffer with 0.09% (w/v) sodium azide and 2% sucrose.
Concentration	0.5 mg/ml
Storage	For short term use, store at 2-8°C up to 1 week. For long term storage, store at -2°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	Synthetic peptide located within the following region: VPVATNDPDEEVQEETRGRFHLLWDPRRKNCLSIRDARRRDNAAYFFRL
Clonality	Polyclonal
MW	47kDa
Uniprot ID	O43699
NCBI	NP_001236
Expiration Date	12 months from date of receipt.

SIGLEC6



Western blot analysis of human Fetal Mus...



Western blot analysis of Jurkat cell lys...